

AMRAD MEDICAL



AC400 | AmRad™ Classic


General Radiographic Systems

www.AmRadMedical.com

Features & Benefits



Perform all exam types to include cross table laterals



Compatible with most DR Solutions

Floor to Wall/Ceiling Tubestand

- Economical yet full-featured, the AC400 can be configured to meet the needs of any department
- Cross-table exams are easily performed using the +/- 180° rotational and +/- 5 in. transverse positioning capabilities
- Available with either 8 ft. or optional 10 ft. rail, the AC400 can be used to perform most weight bearing studies

- ◆ **Compact design for facilities with space constraints**
- ◆ **Economical and versatile to meet the budgetary requirements of the facility**
- ◆ **Customized and flexible configurations**
- ◆ **Industry Leading Warranty**



Completely Flat Table Top



Lowers to Perform Standing Knees

Four-way Float Top Table

- Flat, four-way, floating top permits easy transfer and positioning of patients
- Designed for all types of applications, available in 400 lb patient capacity
- Recessed foot-treadle allows for safe patient transfer

Wallstand

- Designed for ease of use and simplicity, the wallstand lowers for weight bearing exams like standing knees
- Ergonomically designed handle allows for easy positioning of the wall receptor
- The wallstand can be used with CR and DR panels


The Urgent Care Solution



Elevating Table, 400 lb Patient Capacity



Elevating Table, 650 lb Patient Capacity



J1000 Wallstand with PA Patient Hand Grips



Overhead Grip for the J1000 Wallstand

Customize the AC400 System with Options and Accessories (L to R, above)

The AmRad™ Medical AC400 System can be configured to meet the needs of any department. Options include, but are not limited to:

- Elevating Table, 400 lb Patient Capacity
- Elevating Table, 650 lb Patient Capacity
- J1000 Wallstand with PA Hand Grips
- Overhead Grip for J1000 - When Your Patient Needs More Stability
- Grid cabinets in place of Buckys
- 10 ft. Floor Rail tubestand



7555 N. Caldwell Avenue
Niles, Illinois 60714

Phone: 773-353-4000
Fax: 773-588-3424
Toll Free: 1-800-729-9729

AMRAD  **MEDICAL**

www.AmRadMedical.com

find us on 



SAI GLOBAL
ISO 13485
Quality