


AMRAD MEDICAL



DFMT | FMT

Premier Radiographic Systems

www.AmRadMedical.com

* Overhead Grip
on Wallstand, optional




DFMT

Fast and Accurate Patient Positioning

Digital Floor Mounted Tubestand

- Easy to read LCD shows SID continuous to table-top and table receptor with discreet SID of 40" and 72" to the wallstand
- Cross-table views are easily accomplished using the rotational and transverse positioning capabilities
- 10 ft. floor track allows for off-table and weight bearing exams and is customizable down to 8 ft.



Compatible with most DR Solutions



Easy Patient Transfer

Convenient Lock-Release Pedals



Standing Knees



Ergonomic Design

Elevating Four-way Float Top Table

- Flat, four-way, floating top permits easy transfer and positioning of patients
- Designed for all types of applications, available in 400lb and 650lb patient capacities
- Four-pedal design allows for easy accessibility to table lock release
- Collision-protection electronics and fail-safe locks ensure safe operation

Options include:

- 35.5" wide table top
- Four-way lock release handle
- Non-elevating version available in 400lb and 650lb patient capacities

Wallstand

- Center mounting or side mounting of Bucky makes this wallstand flexible for any patient application or installation limitation
- Ergonomically designed hand control for ease of positioning wall receptor
- Side-mounted handgrips stabilize patients during PA exams

Options include:

- Overhead Grip
- Grid Cabinets in place of Bucky



35.5" Wide Table Top



Economical Alternative FMT



Table Lock Release Handle with Elevation



Overhead Grip for J1000 Wallstand

System Options (L to R, above)

The AmRad™ Medical FMT and DFMT Systems can be configured to meet the needs of any department. Options include, but are not limited to:

- An extra-wide table top for the S222 elevating table to help patients feel more secure
- Four-way lock release handle for easy patient positioning can be located anywhere on the rails of the S223 non-elevating or S222 elevating tables
- The FMT is an economical alternative to the DFMT with the same versatility and functionality as the DFMT
- Patient Overhead Grip for use with the J1000 wallstand helps steady patients for upright exams



7555 N. Caldwell Avenue
Niles, Illinois 60714

Phone: 773-353-4000
Fax: 773-588-3424
Toll Free: 1-800-729-9729



www.AmRadMedical.com



SAI GLOBAL
ISO 13485
Quality